Immunoprecipitations

Protocol:

- 1. Use 300-400 μ g protein per IP
- 2. Add primary antibody as directed (usually $\sim 2 \mu g$ for a monoclonal antibody)
- 3. Incubate 1 to 2 hr on ice
- 4. For polyclonals: Add 25 μ l Pansorbin (10% suspension), washed and preincubated in RIPA buffer containing 10 mg/ml BSA
- 5. For monoclonals: Wash and preincubate Pansorbin (10% suspension) in RIPA buffer containing 10 mg/ml BSA. Add 1 μl rabbit anti-mouse IgG/25 μl Pansorbin (10% suspension) and incubate on ice for 45 min. Spin down Pansorbin and remove supernatnat. Wash once in IP buffer containing 10 mg/ml BSA. Make up to original suspension volume in RIPA buffer and add to IP's.
- 6. Incubate IP's with Pansorbin for 30 min on ice.
- 7. Spin down in a microfuge for 10 to 15 sec and aspirate the supernatant
- 8. Wash pellet three times by resuspending in RIPA buffer, pelleting, and removing the supernatant.
- 9. Resuspend final pellet in 50 μ l RIPA buffer and add 50 μ l SDS sample buffer.
- 10. Boil 3 min, pellet the Pansorbin by centrifugation in a microfuge for 30 sec
- 11. Load the supernatant onto an SDS gel

Note:

You may substitute anti-mouse IgG-agarose or protein G-agarose for the Pansorbin plus rabbit-anti-mouse antibodies for use with monoclonal primary antibodies. Use according to manufacturer's instructions.