

Precipitation Lowry Protein Assay (Peterson, Anal. Biochem. 83, 346 (1977))

Solutions Needed:

0.15% Deoxycholate in dH₂O
72% Trichloroacetic acid in dH₂O
20% Na₂CO₃
0.2% CuSO₄ (pentahydrate) + 0.4% KTartrate
0.8 N NaOH
10% SDS
Folin reagent

Make up fresh each time:

CTC Reagent:

10% Na₂CO₃
0.1% CuSO₄
0.2% KTartrate

This solution can be made by diluting equal parts of the Na₂CO₃ stock solution and the 0.2% CuSO₄ (pentahydrate) + 0.4% KTartrate solution.

Lowry Mix:

Add equal volumes of
CTC reagent
0.8 N NaOH
10% SDS
dH₂O

Diluted Folin Reagent:

1 part Folin reagent plus 5 parts dH₂O

Procedure:

1. Aliquot protein samples into 12 x 75 mm plastic tubes
2. To each tube add water to make up to 1 ml volume
3. Add 100 μ l 0.15% DOC
4. Add 100 μ l 72% TCA to precipitate protein
5. Pellet the precipitate in a table top centrifuge by spinning for 30 min at 3000 rpm
6. Pour off supernatant and aspirate any remaining droplets of liquid
7. To each tube add:
 - 400 μ l Lowry mix
 - 400 μ l dH₂O
 - 200 μ l diluted Folin reagent
8. Incubate 30 min at room temperature to allow color development
9. Read samples at 750 nm